

# TM40L rabbit pAb

Catalog No: #44049



Package Size: #44049-1 50ul #44049-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	TM40L rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB
Species Reactivity	Human, Mouse,Rat
Immunogen Description	Synthesized peptide derived from human TM40L
Accession No.	Swiss-Prot#: Q969M1
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

WB 1αO 500-2000

## Images



Western blot analysis of lysates from HeLa cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only and is not intended for use in humans or animals.