

Mouse Small nuclear ribonucleoprotein antibody (snRNP/Sm-Ab) ELISA Kit

Catalog No: #EK20357

Package Size: #EK20357-1 48T #EK20357-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Small nuclear ribonucleoprotein antibody (snRNP/Sm-Ab) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:6.25-400 mU/mL

Sensitivity:2.57mU/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.