

Recombinant Zika virus Genome polyprotein

Catalog No: #AP89643



Package Size: #AP89643-1 20ug #AP89643-2 100ug #AP89643-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Zika virus Genome polyprotein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 85% as determined by SDS-PAGE.
Target Name	His-tagged
Accession No.	Q32ZE1
Target Species	Zika virus (ZIKV)
SDS-PAGE MW	28.2 kDa
Tag Info	N-terminal 6xHis-tagged
Target Sequence	Partial 1415-1463aa+GGGGSGGGG+1499-1668aa KSVDMYIERAGDITWEKDAEVTGNSPRLDVALDESGDFSLVEEDGPPMRGGGGSGGGGSGALWDVPAPKE VKKGETTDGVYRVMTRRLGSGTQVGVGMQEGVFHTMWHVTKGAALRSGEGRDPYWGDKQDLVSYCG PWKLDAAWDGLSEVQLLAVPPGERARNIQLPGIFKTKDGDIGAVALDYPAGTSGSPILDKCGRVIGLYGNGV VIKNGSYVSAITQGKR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.