

Rabbit Alkaline Phosphatase, Liver/Bone/Kidney (ALPL) ELISA Kit

Catalog No: #EK20358



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #EK20358-1 48T #EK20358-2 96T

Description

Product Name	Rabbit Alkaline Phosphatase, Liver/Bone/Kidney (ALPL) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Rabbit

Application Details

Detect Range:0.31-20 ng/mL

Sensitivity:0.135ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.