

Recombinant Mouse Protein disulfide-isomerase A3(Pdia3),Yeast

Catalog No: #AP89644

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP89644-1 20ug #AP89644-2 100ug #AP89644-3 1mg

Description

Product Name	Recombinant Mouse Protein disulfide-isomerase A3(Pdia3),Yeast
Brief Description	Recombinant
Host Species	Yeast
Purification	>85% (SDS-PAGE)
Immunogen Description	Expression Region: 25-505aa Sequence Info:Full Length of Mature Protein
Target Name	His-tagged
Other Names	Pdia3; Erp; Erp60; Grp58; Protein disulfide-isomerase A3; EC 5.3.4.1; 58 kDa glucose-regulated protein; 58 kDa microsomal protein; p58; Disulfide isomerase ER-60; Endoplasmic reticulum resident protein 57; ER protein 57; ERp57; Endoplasmic reticulum resident protein 60; ER protein 60; ERp60
Accession No.	P27773
Target Species	Mus musculus (Mouse)
Tag Info	N-terminal His-tagged
Target Sequence	SDVLEL TDENFESRVS DTGSAGLMLV EFFAPWCGHC KRLAPEYEAA ATRLKGIVPL AKVDCTANTN TCNKYGVSGY PTLKIFRDGE EAGAYDGPRT ADGIVSHLKK QAGPASVPLR TEEEFKFKFIS DKDASVVGFF RDLFSDGHSE FLKAASNLRD NYRFAHTNIE SLVKEYDDNG EGITIFRPLH LANKFEDKTV AYTEKKMTSG KIKKFIQDSI FGLCPHMTED NKDLIQGKDL LTAYYDVDYE KNAKGSNYWR NRVMVAKKF LDAGHKLNFA VASRKTFSHE LSDFGLESTT GEVPVVAIRT AKGEKFVMQE EFSRDGKALE QFLQEYFDGN LKRYLKSEPI PESNEGPVKV VVAENFDDIV NEEDKDVLIIE FYAPWCGHCK NLEPKYKELG EKLSKDPNIV IAKMDATAND VPSPYEVKGF PTIYFSPANK KLTPKKYEGG RELNDFISYL QREATNPPII QEEKPKKKKK AQEDL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.