

Mouse Interferon-induced helicase C domain-containing protein 1 (IFIH1) ELISA Kit

Catalog No: #EK20365

Package Size: #EK20365-1 48T #EK20365-2 96T



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Mouse Interferon-induced helicase C domain-containing protein 1 (IFIH1) ELISA Kit |
| Brief Description | ELISA Kit |
| Applications | ELISA |
| Species Reactivity | Mouse |

Application Details

Detect Range:15.6-1000pg/mL

Sensitivity:6.2pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.