

# Recombinant Human Phosphatidylinositol 3-kinase catalytic subunit type 3(PIK3C3) ,partial

Catalog No: #AP89645

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP89645-1 20ug #AP89645-2 100ug #AP89645-3 1mg

## Description

Product Name	Recombinant Human Phosphatidylinositol 3-kinase catalytic subunit type 3(PIK3C3) ,partial
Brief Description	Recombinant
Host Species	E.coli
Purification	>85% (SDS-PAGE)
Immunogen Description	Sequence Info:Partial
Target Name	His-tagged
Other Names	hVps34; MGC61518; Phosphatidylinositol 3 kinase catalytic subunit type 3; Phosphatidylinositol 3 kinase class 3; Phosphatidylinositol 3 kinase p100 subunit; Phosphatidylinositol 3-kinase catalytic subunit type 3; Phosphatidylinositol 3-kinase p100 subunit; Phosphoinositide 3 kinase class 3; Phosphoinositide-3-kinase class 3; PI3 kinase type 3; PI3-kinase type 3; PI3K type 3; Pik3c3; PK3C3_HUMAN; PtdIns 3 kinase type 3; PtdIns-3-kinase type 3; Vps 34; Vps34
Accession No.	Q8NEB9
Target Species	Homo sapiens (Human)
Tag Info	N-terminal His-tagged
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.