

Rat Ras Homolog Gene Family, Member A (RHOA) ELISA Kit

Catalog No: #EK20367

Package Size: #EK20367-1 48T #EK20367-2 96T



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Ras Homolog Gene Family, Member A (RHOA) ELISA Kit
Brief Description	Rat Ras Homolog Gene Family, Member A (RHOA) ELISA Kit
Applications	ELISA
Species Reactivity	Rat

Application Details

Detect Range:0.156-10ng/mL

Sensitivity:0.055ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.