

## Recombinant Human Plasminogen(PLG)

Catalog No: #AP89646

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Plasminogen(PLG)
Brief Description	Recombinant
Host Species	E.coli
Purification	>85% (SDS-PAGE)
Immunogen Description	Expression Region: 20-810aa Sequence Info:Full Length of Mature Protein
Target Name	His-tagged
Other Names	Plasmin; Plasmin heavy chain A; Plasmin light chain B; Plasminogen; PLG; PLMN_HUMAN
Accession No.	P00747
Target Species	Homo sapiens (Human)
Tag Info	N-terminal His-tagged
Target Sequence	E PLDDYVNTQG ASLFSVTKKQ LGAGSIEECA AKCEEDEF CRAFTQYHSKE QQCIVMAENR KSSIIIRMRD VVLFEKKVYL SECKTGNNGKN YRGTMSKTKN GITCQKWSST SPHRPRFSPA THPSEGLEEN YCRNPNDPQ GPWCYTTDPE KRYDYCDILE CEEECMHCSG ENYDGKISK TMSGLECAWD SQSPHAHGYI PSKFPNKNL KNYCRNPDR LRPWCFTTDP NKRWELCDIP RCTTPPPSSG PTYQCLKGTG ENYRGNVAVT VSGHTCQHWS AQTPTHNRT PENFPCKNLD ENYCRNPDGK RAPWCHTNS QVRWEYCKIP SCDSSPVSTE QLAPTAPPEL TPVVQDCYHG DGQSYRGTS TTTTGKKCQS WSSMTPHRHQ KTPENYPNAG LTMNYCRNP ADKGPWCFTT DPSVRWEYCN LKKCSGTEAS VVAPPPVLL PDVETPSEED CMFGNGKGYR GKRATTVTGT PCQDWAQEP HRHSIFTPET NPRAGLEKNY CRNPDGDVGG PWCYTTNPRK LYDYCDVPQC AAPSFDCGKP QVEPKKCPGR VGGCVAHPH SWPWQVSLRT RFGMHFCGGT LISPEWVLT AHCLEKSPRP SSYKVLGAH QEVNLEPHVQ EIEVSRLFLE PTRKDIALK LSSPAVITDK VIPACLPSPN YVVADRTECF ITGWGETQGT FGAGLLKEAQ LPVIENKVCN RYEFNLGRVQ STELCAGHLA GGTDCSQGDS GGPLVCFEKD KYILQGVTSW GLGCARPKNP GYVVRVSRFV TWIEGVMRNN
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.