

# Mouse Glycoprotein Ib Alpha Polypeptide, Platelet (GP1Ba) ELISA Kit

Catalog No: #EK20376



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #EK20376-1 48T #EK20376-2 96T

## Description

Product Name	Mouse Glycoprotein Ib Alpha Polypeptide, Platelet (GP1Ba) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

## Application Details

Detect Range: 0.156-10 ng/mL

Sensitivity: 0.053 ng/mL

Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50  $\mu$ L

Assay Time: 1-4.5h

Detection wavelength: 450 nm

## Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.