

TREM-1 Polyclonal Antibody

Catalog No: #44120



Package Size: #44120-1 50ul #44120-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	TREM-1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IHC-p;IF/ICC
Species Reactivity	Human;Mouse;Rat
Specificity	TREM-1 Polyclonal Antibody detects endogenous levels of TREM-1 protein.
Immunogen Description	Synthesized peptide derived from TREM-1 .
Accession No.	Swiss-Prot#: Q9NP99 NCBI Gene ID:54210
Calculated MW	26KD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1:500-2000

IHC-p 1:100-500

Images



Western Blot analysis of THP-1, mouse heart cells using TREM-1 Polyclonal Antibody.



Immunohistochemical analysis of paraffin-embedded human-brain, antibody was diluted at 1:100

Note: This product is for in vitro research use only and is not intended for use in humans or animals.