

Human Pre-mRNA-processing factor 40 homolog A (PRPF40A) ELISA Kit

Catalog No: #EK20389

Package Size: #EK20389-1 48T #EK20389-2 96T



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human Pre-mRNA-processing factor 40 homolog A (PRPF40A) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:78.1-5000pg/mL
Sensitivity:35.3pg/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.