

Recombinant Mouse Protachykinin-1(Tac1)

Catalog No: #AP89649



Package Size: #AP89649-1 20ug #AP89649-2 100ug #AP89649-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse Protachykinin-1(Tac1)
Brief Description	Recombinant Proteins
Host Species	Yeast
Purification	>85% (SDS-PAGE)
Immunogen Description	Expression Region: 58-130aa
Target Name	His-tagged
Other Names	Tac1; Nka; Nkna; Tac2; Protachykinin-1; PPT) [Cleaved into: Substance P; Neurokinin A; NKA; Neuromedin L; Substance K); Neuropeptide K; NPK); Neuropeptide gamma; C-terminal-flanking peptide]
Accession No.	P41539
Target Species	Mus musculus (Mouse)
Tag Info	N-terminal His-tagged
Target Sequence	Sequence Info:Full Length of Mature Protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.