

POGK rabbit pAb

Catalog No: #44136



Package Size: #44136-1 50ul #44136-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	POGK rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB;IHC-p
Species Reactivity	Human; Mouse
Specificity	This antibody detects endogenous levels of POGK at Human/Mouse
Immunogen Description	Synthesized peptide derived from human POGK
Accession No.	Swiss Prot:Q9P215GeneID:57645
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


WB 1:500-2000

IHC-p 1:50-300

Images



Western blot analysis of lysates from 3T3 cells, primary antibody was diluted at 1:1000, 4° over night



Immunohistochemical analysis of paraffin-embedded human tonsil.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.