

Rabbit Syndecan-1 (SDC1) ELISA Kit

Catalog No: #EK20421



Package Size: #EK20421-1 48T #EK20421-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rabbit Syndecan-1 (SDC1) ELISA Kit
Host Species	Rabbit
Applications	ELISA

Application Details

Detect Range:31.25-2000 pg/mL

Sensitivity:0.66ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:competitive inhibition

Note: This product is for in vitro research use only and is not intended for use in humans or animals.