

# Streptavidin AF488 Conjugated

Catalog No: #LS1013



Package Size: #LS1013-1 100ul #LS1013-2 1ml

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Streptavidin AF488 Conjugated
Purification	>90% as determined by SDS-PAGE
Conjugates	AF488
Concentration	oO 0.5 mg/ml
Formulation	10 mM TBS (pH=7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.