

# Rabbit Colony Stimulating Factor 3, Granulocyte (GCSF) ELISA Kit

Catalog No: #EK20435



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #EK20435-1 48T #EK20435-2 96T

## Description

Product Name	Rabbit Colony Stimulating Factor 3, Granulocyte (GCSF) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Rabbit

## Application Details

Detect Range:7.8-500pg/mL

Sensitivity:3.1pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.