

Rabbit Colony Stimulating Factor 2, Granulocyte Macrophage (GMCSF) ELISA Kit

Catalog No: #EK20452



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #EK20452-1 48T #EK20452-2 96T

Description

Product Name	Rabbit Colony Stimulating Factor 2, Granulocyte Macrophage (GMCSF) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Rabbit

Application Details

Detect Range:15.6-1000pg/mL

Sensitivity:7.2pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.