

Human 28S ribosomal protein S34, mitochondrial (MRPS34) ELISA Kit

Catalog No: #EK20441



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #EK20441-1 48T #EK20441-2 96T

Description

Product Name	Human 28S ribosomal protein S34, mitochondrial (MRPS34) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:125-8000pg/mL
Sensitivity:58.8pg/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.