

## MUC2 (ABT-MUC2) mouse mAb

Catalog No: #44168

Package Size: #44168-1 50ul #44168-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	MUC2 (ABT-MUC2) mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Applications	IHC-p;IF
Species Reactivity	Human
Specificity	This antibody detects endogenous levels of human MUC2. Heat-induced epitope retrieval (HIER) TRIS-EDTA of pH8.0 was highly recommended as antigen repair method in paraffin section
Immunogen Description	Synthesized peptide derived from human MUC2
Accession No.	Swiss-Prot#: Q02817; NCBI Gene#: 4583
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

IHC-p 1:100-500

## Images



Immunohistochemical analysis of paraffin-embedded Colon.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.