

Human CDH15/ Cadherin-15 ELISA Kit

Catalog No: #EK20445



Package Size: #EK20445-1 48T #EK20445-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human CDH15/ Cadherin-15 ELISA Kit
Brief Description	Human CDH15/ Cadherin-15 ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.312-20ng/ml

Sensitivity:0.125ng/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.