

Rat Runt Related Transcription Factor 2 (RUNX2) ELISA Kit

Catalog No: #EK20456

Package Size: #EK20456-1 48T #EK20456-2 96T



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Runt Related Transcription Factor 2 (RUNX2) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Rat
Other Names	RUNX2

Application Details

Detect Range:0.312-29ng/ml

Sensitivity:0.071ng/ml

Sample Type:serum, plasma or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.