

Erk1/2(Phospho-T202/Y204+T185/Y187) Rabbit mAb

Catalog No: #14151



Package Size: #14151-1 50ul #14151-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Erk1/2(Phospho-T202/Y204+T185/Y187) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human Mouse Rat
Specificity	Phospho-ERK1 (Y204) + ERK2 (Y187) Antibody detects endogenous levels of total Phospho-ERK1 (Y204) + ERK2 (Y187)
Immunogen Description	A synthesized peptide derived from human Phospho-ERK1 (Y204) + ERK2 (Y187)
Other Names	ERK-1; ERK1; ERT2; kinase ERK1; MAP kinase 1; MAPK 1; MAPK3; MK03; MNK1; p44-ERK1; P44-ERK1; p44-MAPK;
Accession No.	Uniprot:P27361/P28482
Calculated MW	42kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IP:1:50

Images



Western blot analysis of Phospho-Erk1 (T202/Y204) + Erk2 (T185/Y187) expression in A431 cell lysate treated with EGF.

Product Description

Serine/threonine kinase which acts as an essential component of the MAP kinase signal transduction pathway. MAPK1/ERK2 and MAPK3/ERK1 are the 2 MAPKs which play an important role in the MAPK/ERK cascade. They participate also in a signaling cascade initiated by activated KIT and KITLG/SCF. Depending on the cellular context, the MAPK/ERK cascade mediates diverse biological functions such as cell growth, adhesion, survival

and differentiation through the regulation of transcription, translation, cytoskeletal rearrangements.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.