

JNK1/2/3 (Phospho-Thr183+Thr183+Thr221) Rabbit mAb

Catalog No: #14169

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #14169-1 50ul #14169-2 100ul


Description

Product Name	JNK1/2/3 (Phospho-Thr183+Thr183+Thr221) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human Mouse Rat
Specificity	Phospho-JNK1/2/3 (T183+T183+T221) Antibody detects endogenous levels of total Phospho-JNK1/2/3 (T183+T183+T221)
Immunogen Description	A synthesized peptide derived from human Phospho-JNK1/2/3 (T183+T183+T221)
Other Names	JNK 46; JNK 55; MAPK10; MAPK9; MAPK8; SAPK1b ; SAPK1; SAPK; PRKM10; PRKM9; PRKM8;
Accession No.	Uniprot:P45983/P45984/P53779
Calculated MW	46,54kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200IP:1:50


Images



Western blot analysis of JNK1/2/3 phosphorylation expression in NIH/3T3 cell lysate treated with Anisomycin.



Immunohistochemical analysis of paraffin-embedded human breast, using Phospho-JNK1/2/3 (T183+T183+T221) Antibody.



Immunofluorescent analysis of NIH/3T3 cells treated with Anisomycin, using Phospho-JNK1/2/3 (T183+T183+T221) Antibody.

Product Description

Phospho-JNK1/2/3 (T183+T183+T221) Antibody will detect will detect JNK1 (pT183), JNK2 (pT183) and JNK3 (pT221).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.