

## c-Jun (Phospho-Thr91) Rabbit mAb

Catalog No: #14214

Package Size: #14214-1 50ul #14214-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	c-Jun (Phospho-Thr91) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human Mouse Rat
Specificity	Phospho-c-Jun (T91) Antibody detects endogenous levels of Phospho-c-Jun (T91)
Immunogen Description	A synthesized peptide derived from human c-Jun
Other Names	JUN; p39; V-jun avian sarcoma virus 17 oncogene homolog; Proto-oncogene c-Jun; AP1; Activator protein 1; Transcription factor AP-1;
Accession No.	Uniprot:P05412
Calculated MW	43kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:1000~1:2000IP:1:20

## Images



Western blot analysis of Phospho-c-Jun (T91) expression in human brain lysate.

## Product Description

c-Jun is a member of the Jun family containing c-Jun, JunB, and JunD, and is a component of the transcription factor activator protein-1 (AP-1). AP-1 regulated genes exert diverse biological functions including cell proliferation, differentiation, and apoptosis, as well as transformation, invasion and metastasis, depending on cell type and context.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.