

PAK4 + PAK5 + PAK6 (S474 + S560 + S602) Rabbit mAb

Catalog No: #14252

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #14252-1 50ul #14252-2 100ul


Description

Product Name	PAK4 + PAK5 + PAK6 (S474 + S560 + S602) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human Mouse Rat
Specificity	Phospho-PAK4 + PAK5 + PAK6 (S474 + S560 + S602) Antibody detects endogenous levels of total Phospho-PAK4 + PAK5 + PAK6 (S474 + S560 + S602)
Immunogen Description	A synthesized peptide derived from human Phospho-PAK4 + PAK5 + PAK6 (S474 + S560 + S602)
Other Names	PAK4; PAK5; PAK6;
Accession No.	Uniprot:O96013(PAK4-Antigen),Q9P286(PAK5),Q9NQU5(PAK6)
Calculated MW	64, 80, 75 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000FC:1:200

Images



Western blot analysis of Phospho-PAK4 + PAK5 + PAK6 (S474 + S560 + S602) expression in U87-MG cell lysate.

Product Description

Serine/threonine protein kinase that plays a role in a variety of different signaling pathways including cytoskeleton regulation, cell migration, growth, proliferation or cell survival.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.