

# Bad (Phospho-Ser112) Rabbit mAb

Catalog No: #14254

Package Size: #14254-1 50ul #14254-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	Bad (Phospho-Ser112) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human Mouse Rat
Specificity	Phospho-Bad (S112) Antibody detects endogenous levels of total Phospho-Bad (S112)
Immunogen Description	A synthesized peptide derived from human Bad
Other Names	BAD; BBC2; BBC6; BCL2L8;
Accession No.	Uniprot:Q92934
Calculated MW	23kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000IP:1:20

## Images



Western blot analysis of Phospho-Bad (S112) expression in (1) HeLa cell lysate; (2) HeLa cell treated with Calcyculin A lysate.

## Product Description

Promotes cell death. Successfully competes for the binding to Bcl-X(L), Bcl-2 and Bcl-W, thereby affecting the level of heterodimerization of these proteins with BAX.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.