

TAK1(Phospho-Ser439) Rabbit mAb

Catalog No: #14266



Package Size: #14266-1 50ul #14266-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TAK1(Phospho-Ser439) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human Mouse Rat
Specificity	Phospho-TAK1 (S439) Antibody detects endogenous levels of total Phospho-TAK1 (S439)
Immunogen Description	A synthesized peptide derived from human Phospho-TAK1 (S439)
Other Names	MAP3K 7; MEKK7; Mitogen activated protein kinase kinase kinase 7; TAK1; TGF beta activated kinase 1; TGF1a;
Accession No.	Uniprot:O43318
Calculated MW	72kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IP:1:50

Images



Western blot analysis of Phospho-TAK1 (S439) expression in (1) HeLa treated with Calyculin A + IL-1 beta cell lysate; (2) Untreated.

Product Description

Component of a protein kinase signal transduction cascade. Mediator of TRAF6 and TGF-beta signal transduction. Activates IKBKB and MAPK8 in response to TRAF6 signaling. Stimulates NF-kappa-B activation and the p38 MAPK pathway. In osmotic stress signaling, plays a major role in the activation of MAPK8/JNK, but not that of NF-kappa-B.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.