

## Myosin heavy chain Rabbit mAb

Catalog No: #56130



Package Size: #56130-1 50ul #56130-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	Myosin heavy chain Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse Rat
Specificity	Myosin heavy chain Antibody detects endogenous levels of total Myosin heavy chain
Immunogen Description	A synthesized peptide derived from human Myosin heavy chain
Other Names	MYH6; MYH7; MYHC; MYHCA; MYHCB; Myosin heavy chain 6; Myosin heavy chain; Myosin-6;
Accession No.	Uniprot:P12883/P13533
Calculated MW	224kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


## Application Details

WB:1:500~1:1000IHC:1:50~1:200

## Images



Western blot analysis of Myosin heavy chain expression in mouse heart lysate.



Immunohistochemical analysis of paraffin-embedded human skeletal muscle, using Myosin heavy chain Antibody.

## Product Description

---

Muscle contraction.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.