

## Substance P Rabbit mAb

Catalog No: #56143



Package Size: #56143-1 50ul #56143-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	Substance P Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	IHC
Species Reactivity	Human
Specificity	Substance P Antibody detects endogenous levels of total Substance P
Immunogen Description	A synthesized peptide derived from human Substance P
Other Names	Protachykinin-1; PPT; Neurokinin A; Neuromedin L; Substance K; Neuropeptide K; NPK; TAC1; NKA, NKNA, TAC2; Tachykinin 1; tachykinin 2;
Accession No.	Uniprot:P20366
Calculated MW	kDa
Storage	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

IHC:1:50~1:200

## Images



Immunohistochemical analysis of paraffin-embedded human brain, using Substance P Antibody.

## Product Description

Tachykinins are active peptides which excite neurons, evoke behavioral responses, are potent vasodilators and secretagogues, and contract (directly or indirectly) many smooth muscles.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.