

## FGFR4 Rabbit mAb

Catalog No: #56179

Package Size: #56179-1 50ul #56179-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	FGFR4 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	FGFR4 Antibody detects endogenous levels of total FGFR4
Immunogen Description	A synthesized peptide derived from human FGFR4
Other Names	Fibroblast growth factor receptor 4; FGFR-4; CD334; FGFR4; JTK2; TKF;
Accession No.	Uniprot:P22455
Calculated MW	88kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB: 1:500~1:2000

## Images



Western blot analysis of FGFR4 expression in Raji cell lysate.

## Product Description

Receptor for acidic fibroblast growth factor. Does not bind to basic fibroblast growth factor. Binds FGF19.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.