

PAI1 Rabbit mAb

Catalog No: #56181

Package Size: #56181-1 50ul #56181-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PAI1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP FC
Species Reactivity	Human
Specificity	PAI1 Antibody detects endogenous levels of total PAI1
Immunogen Description	A synthesized peptide derived from human PAI1
Other Names	PAI-1; PAI-1; PAI; Endothelial plasminogen activator inhibitor; SERPINE1; PAI1; PLANH1; Serpine 1;
Accession No.	Uniprot:P05121
Calculated MW	45kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:1000 ICC/IF:1:50~1:200 IP:1:30 FC:1:50

Images



Western blot analysis of PAI1 expression in HepG2 cell lysate.

Product Description

This inhibitor acts as 'bait' for tissue plasminogen activator, urokinase, and protein C. Its rapid interaction with TPA may function as a major control point in the regulation of fibrinolysis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.