

## FABP3 Rabbit mAb

Catalog No: #56217

Package Size: #56217-1 50ul #56217-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	FABP3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse Rat
Specificity	FABP3 Antibody detects endogenous levels of total FABP3
Immunogen Description	A synthesized peptide derived from human FABP3
Other Names	FABP 11; FABP 3; H FABP; heart; M FABP; MDGI;
Accession No.	Uniprot:P05413
Calculated MW	15kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000IHC:1:50~1:200

## Images



Western blot analysis of FABP3 expression in human fetal muscle lysate.

## Product Description

FABP are thought to play a role in the intracellular transport of long-chain fatty acids and their acyl-CoA esters.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.