

Parathyroid Hormone Rabbit mAb

Catalog No: #56230



Package Size: #56230-1 50ul #56230-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Parathyroid Hormone Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	Parathyroid Hormone Antibody detects endogenous levels of total Parathyroid Hormone
Immunogen Description	A synthesized peptide derived from human Parathyroid Hormone
Other Names	hPTH; Parathormone; Parathyrin; Parathyroid hormone 1; Prepro PTH; Preproparathyroid hormone; PTH1; PTH1 receptor; PTH1R; PTHR1;
Accession No.	Uniprot:P01270
Calculated MW	36kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000 IHC: 1:50~1:200

Images



Western blot analysis of Parathyroid Hormone expression in human Parathyroid Hormone recombinant protein.

Product Description

PTH elevates calcium level by dissolving the salts in bone and preventing their renal excretion. Stimulates [1-14C]-2-deoxy-D-glucose (2DG) transport and glycogen synthesis in osteoblastic cells.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.