

Chromogranin B Rabbit mAb

Catalog No: #56232

Package Size: #56232-1 50ul #56232-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Chromogranin B Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	Chromogranin B Antibody detects endogenous levels of total Chromogranin B
Immunogen Description	A synthesized peptide derived from human Chromogranin B
Other Names	CCB peptide; CgB; Chgb; ChromograninB; SCG1; Secretogranin 1; Secretogranin B; Sgl;
Accession No.	Uniprot:P05060
Calculated MW	125kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000

Images



Western blot analysis of Chromogranin B expression in (1) HeLa cell lysate; (2) Raw 264.7 cell lysate.

Product Description

Secretogranin-1 is a neuroendocrine secretory granule protein, which may be the precursor for other biologically active peptides.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.