

Trefoil Factor 3 Rabbit mAb

Catalog No: #56250



Package Size: #56250-1 50ul #56250-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

| | |
|-----------------------|-----------------------------------------------------------------------------------------------------------------|
| Product Name | Trefoil Factor 3 Rabbit mAb |
| Host Species | Rabbit |
| Clonality | Monoclonal |
| Isotype | Rabbit IgG |
| Purification | Affinity-chromatography |
| Applications | WB IHC IP |
| Species Reactivity | Human |
| Specificity | Trefoil Factor 3 Antibody detects endogenous levels of total Trefoil Factor 3 |
| Immunogen Description | A synthesized peptide derived from human Trefoil Factor 3 |
| Other Names | hITF; hP1.B; Intestinal trefoil factor; ITF; mITF; P1B; Polypeptide P1.B; TFF3; TFI; TREFOIL; Trefoil factor 3; |
| Accession No. | Uniprot:Q07654 |
| Calculated MW | 9kDa |
| Formulation | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage | Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle. |

Application Details

WB:1:500~1:1000IHC:1:100~1:500IP:1:30

Images



Western blot analysis of Trefoil Factor 3 expression in human small intestine lysate.

Product Description

Involved in the maintenance and repair of the intestinal mucosa. Promotes the mobility of epithelial cells in healing processes (motogen).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.