

RNF8 Rabbit mAb

Catalog No: #56279

Package Size: #56279-1 50ul #56279-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RNF8 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Rat
Specificity	RNF8 Antibody detects endogenous levels of total RNF8
Immunogen Description	A synthesized peptide derived from human RNF8
Other Names	C3HC4 type zinc finger protein; RNF8;
Accession No.	Uniprot:O76064
Calculated MW	60kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000

Images



Western blot analysis of RNF8 expression in human fetal kidney lysate.

Product Description

Mediates the ubiquitination of histones H2A and H2AX, thereby promoting the formation of TP53BP1 and BRCA1 ionizing radiation-induced foci (IRIF). Promotes the formation of 'Lys-63'-linked polyubiquitin chains and functions with the specific ubiquitin-conjugating UBE2N/UBC13.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.