

Factor X Rabbit mAb

Catalog No: #56283

Package Size: #56283-1 50ul #56283-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Factor X Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	Factor X Antibody detects endogenous levels of total Factor X
Immunogen Description	A synthesized peptide derived from human Factor X
Other Names	Activated factor Xa heavy chain; Coagulation factor X; F10; factor Xa; FXA; Prothrombinase; Stuart factor;
Accession No.	Uniprot:P00742
Calculated MW	74kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:1000~1:5000IHC:1:1000~1:500

Images



Western blot analysis of Factor X expression in HepG2 cell lysate.

Product Description

Factor Xa is a vitamin K-dependent glycoprotein that converts prothrombin to thrombin in the presence of factor Va, calcium and phospholipid during blood clotting.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.