

GFP Rabbit mAb

Catalog No: #56315

Package Size: #56315-1 50ul #56315-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GFP Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	All
Specificity	GFP antibody detects endogenous levels of C-terminal, internal, and N-terminal GFP-tagged proteins.
Immunogen Description	A synthesized peptide derived from green fluorescent protein
Other Names	GFP ; yfp ; Green fluorescent protein ;
Accession No.	Uniprot:P42212
Calculated MW	27kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:3000~1:10000

Images



Western blot analysis of GFP expression in 293 cell lysate transfected with GFP.

Product Description

Green fluorescent protein(GFP) has been engineered to produce a vast number of variously colored mutants, fusion proteins, and biosensors. Fluorescent proteins and its mutated allelic forms, blue, cyan and yellow have become a useful and ubiquitous tool for making chimeric proteins, where they function as a fluorescent protein tag. Typically they tolerate N- and C-terminal fusion to a broad variety of proteins.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.