

alpha Tubulin yeast Rabbit mAb

Catalog No: #56316



Package Size: #56316-1 50ul #56316-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	alpha Tubulin yeast Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Saccharomyces cerevisiae yeast
Specificity	alpha Tubulin yeast Antibody detects endogenous levels of total alpha Tubulin yeast
Immunogen Description	A synthesized peptide derived from yeast alpha Tubulin
Other Names	Alpha-tubulin 1;H2 ALPHA;LIS3;TUBA1;TUBA1A;Tuba4a;Tubulin alpha-4A chain
Accession No.	Uniprot:P09733
Calculated MW	50kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:3000~1:10000

Images



Western blot analysis of alpha Tubulin expression in Saccharomyces cerevisiae (yeast) cell lysate.

Product Description

Tubulin is the major constituent of microtubules. It binds two moles of GTP, one at an exchangeable site on the beta chain and one at a non-exchangeable site on the alpha chain.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.