

## VSV-G-tag Rabbit mAb

Catalog No: #56318



Package Size: #56318-1 50ul #56318-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	VSV-G-tag Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	All
Specificity	VSV-G tag antibody detects endogenous levels of C-terminal, internal, and N-terminal VSV-G-tagged proteins.
Immunogen Description	A synthesized peptide
Other Names	VSV epitope tag;VSV tag;VSV-G epitope tag;VSVG;vsvg tag
Accession No.	Uniprot:N/A
Calculated MW	N/A
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:3000~1:10000

## Images



Western blot analysis of extracts from VSV-G tag fusion protein, using VSV-G tag antibody.

## Product Description

Vesicular stomatitis virus (VSV), an enveloped RNA virus from the Rhabdoviridae family, is released from the plasma membrane of host cells by a process called budding. The glycoprotein (VSV-g) contains a domain in its extracellular membrane proximal stem that appears to be needed for efficient VSV budding.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.