

## MUC5B Rabbit mAb

Catalog No: #56320

Package Size: #56320-1 50ul #56320-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	MUC5B Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human
Specificity	MUC5B Antibody detects endogenous levels of total MUC5B
Immunogen Description	A synthesized peptide derived from human MUC5B
Other Names	MG1; MUC5; MUC5B; MUC9; Mucin 5B;
Accession No.	Uniprot:Q9HC84
Calculated MW	596kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000FC:1:100

## Images



Western blot analysis of MUC5B expression in HeLa cell lysate.

## Product Description

Gel-forming mucin that is thought to contribute to the lubricating and viscoelastic properties of whole saliva and cervical mucus.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.