

Netrin 1 Rabbit mAb

Catalog No: #56336

Package Size: #56336-1 50ul #56336-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Netrin 1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	Netrin 1 Antibody detects endogenous levels of total Netrin 1
Immunogen Description	A synthesized peptide derived from human Netrin 1
Other Names	Netrin1; NTN1; NTN1L; Unc6;
Accession No.	Uniprot:O95631
Calculated MW	75kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000 ICC/IF:1:50~1:200

Images



Western blot analysis of Netrin 1 expression in Human fetal brain lysate.

Product Description

Netrins control guidance of CNS commissural axons and peripheral motor axons. Its association with either DCC or some UNC5 receptors will lead to axon attraction or repulsion, respectively.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.