

Retinoic Acid Receptor beta Rabbit mAb

Catalog No: #56343

Package Size: #56343-1 50ul #56343-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Retinoic Acid Receptor beta Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human
Specificity	Retinoic Acid Receptor beta Antibody detects endogenous levels of total Retinoic Acid Receptor beta
Immunogen Description	A synthesized peptide derived from human Retinoic Acid Receptor beta
Other Names	HAP; NR1B2; RAR B; RAR beta; RAR epsilon; RARB; RRB2;
Accession No.	Uniprot:P10826
Calculated MW	55kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200

Images



Western blot analysis of Retinoic Acid Receptor beta expression in MCF7 cell lysate.

Product Description

Receptor for retinoic acid. Retinoic acid receptors bind as heterodimers to their target response elements in response to their ligands, all-trans or 9-cis retinoic acid, and regulate gene expression in various biological processes.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.