

## NKX3.1 Rabbit mAb

Catalog No: #56346

Package Size: #56346-1 50ul #56346-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	NKX3.1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	Nkx3.1 Antibody detects endogenous levels of total Nkx3.1
Immunogen Description	A synthesized peptide derived from human Nkx3.1
Other Names	NKX3; BAPX2; NKX3A; NKX3.1; NKX3-1;
Accession No.	Uniprot:Q99801
Calculated MW	32kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:1000IHC:1:50~1:200

## Images



Western blot analysis of Nkx3.1 expression in LNCaP cell lysate.

## Product Description

Transcription factor, which binds preferentially the consensus sequence 5'-TAAGT[AG]-3' and can behave as a transcriptional repressor. Play an important role in normal prostate development, regulating proliferation of glandular epithelium and in the formation of ducts in prostate.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.