

## IMP3 Rabbit mAb

Catalog No: #56358



Package Size: #56358-1 50ul #56358-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	IMP3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP FC
Species Reactivity	Human Mouse Rat
Specificity	IMP3 Antibody detects endogenous levels of total IMP3
Immunogen Description	A synthesized peptide derived from human IMP3
Other Names	CT98; hKOC; IGF2BP3; IMP 3; KOC1; VICKZ 3;
Accession No.	Uniprot:O00425
Calculated MW	69kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000 ICC/IF:1:50~1:200 IP:1:50 FC:1:50

## Images



Western blot analysis of IMP3 expression in HeLa cell lysate.

## Product Description

RNA-binding protein that act as a regulator of mRNA translation and stability. Binds to the 5'-UTR of the insulin-like growth factor 2 (IGF2) mRNAs. Binds to sequences in the 3'-UTR of CD44 mRNA.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.