

Uroplakin3a Rabbit mAb

Catalog No: #56361



Package Size: #56361-1 50ul #56361-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Uroplakin3a Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	IHC FC
Species Reactivity	Human
Specificity	Uroplakin III Antibody detects endogenous levels of total Uroplakin III
Immunogen Description	A synthesized peptide derived from human Uroplakin III
Other Names	UP3a; UPIII; UPIIIA; UPK3; UPK3A; Uroplakin3A; UroplakinIII;
Accession No.	Uniprot:O75631
Calculated MW	kDa
Storage	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

IHC:1:1000~1:500FC:1:100

Images



Immunohistochemical analysis of paraffin-embedded human bladder, using Uroplakin III Antibody.

Product Description

Component of the asymmetric unit membrane (AUM); a highly specialized biomembrane elaborated by terminally differentiated urothelial cells. May play an important role in AUM-cytoskeleton interaction in terminally differentiated urothelial cells.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.