

RANK Rabbit mAb

Catalog No: #56372

Package Size: #56372-1 50ul #56372-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RANK Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	RANK Antibody detects endogenous levels of total RANK
Immunogen Description	A synthesized peptide derived from human RANK
Other Names	CD265; FEO; LOH18CR1; mRANK; ODFR; OFE; OPTB7; OSTs; PDB2; RANK; TNFRSF11A; TRANCER;
Accession No.	Uniprot:Q9Y6Q6
Calculated MW	80kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000

Images



Western blot analysis of RANK expression in Human fetal liver lysate.

Product Description

Receptor for TNFSF11/RANKL/TRANCE/OPGL; essential for RANKL-mediated osteoclastogenesis. Involved in the regulation of interactions between T-cells and dendritic cells.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.