

Thyroglobulin Rabbit mAb

Catalog No: #56382

Package Size: #56382-1 50ul #56382-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Thyroglobulin Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human
Specificity	Thyroglobulin Antibody detects endogenous levels of total Thyroglobulin
Immunogen Description	A synthesized peptide derived from human Thyroglobulin
Other Names	AITD3; hTG; TDH3; Tg; TGN; Thyroglobulin;
Accession No.	Uniprot:P01266
Calculated MW	305kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:1000~1:5000IHC:1:50~1:200ICC/IF:1:50~1:200IP:1:50

Images



Western blot analysis of Thyroglobulin expression in Human thyroid lysate.

Product Description

Precursor of the iodinated thyroid hormones thyroxine (T4) and triiodothyronine (T3).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.