

PDE4D Rabbit mAb

Catalog No: #56424

Package Size: #56424-1 50ul #56424-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PDE4D Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	PDE4D Antibody detects endogenous levels of total PDE4D
Immunogen Description	A synthesized peptide derived from human PDE4D
Other Names	DPDE3; HSPDE4D; PDE43; Pde4d; PDE4DN2; STRK1;
Accession No.	Uniprot:Q08499
Calculated MW	91kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000

Images



Western blot analysis of PDE4D expression in K562 cell lysate.

Product Description

Hydrolyzes the second messenger cAMP, which is a key regulator of many important physiological processes.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.